

## Counting by Measuring Mass Molar Mass Lab

**Purpose:** To determine the mass of several samples of chemical compounds and use the data to count atoms

**Materials:** H<sub>2</sub>O(l) NaCl(s) CaCO<sub>3</sub>(s) plastic spoon  
Watch glass small beaker balance paper  
Pencil ruler weighing paper

**Procedure:** Measure the mass of one level teaspoon of sodium chloride, water, and calcium carbonate. Place your results in the table.

**Analyze:** Use your data to complete the following steps. Record your answers in the data table.

1. Calculate the moles of NaCl contained in one level teaspoon. Use the periodic table to determine the molar mass of sodium chloride.
2. Repeat Step 1 for the other two compounds. Use the periodic table to calculate the molar mass of water and calcium carbonate.
3. Calculate the number of moles of each element present in a teaspoon sized sample of H<sub>2</sub>O. Use dimensional analysis. Repeat for the other two compounds.
4. Calculate the number of atoms of each element present in the teaspoon sized sample of H<sub>2</sub>O. Use Avogadro's number. Repeat for the other compounds in your table.
5. Which of the three teaspoon sized samples contains the greatest number of moles? \_\_\_\_\_
6. Which of the three compounds contains the most atoms? \_\_\_\_\_

**DATA TABLE:**

	H <sub>2</sub> O(l)	NaCl(s)	CaCO <sub>3</sub> (s)
Mass(grams)			
Molar Mass(g/mol)			
Moles of each compound			
Moles of each element			
Atoms of each element			

**All Calculations(show below):**